How to use your new Rehband Support - 6903

UK

Function

Offers support to help relieve pain in the area.

Indication

Overuse, inflammation in hand/underarm, light knee instabilities

Material

Polyamide, polyester/bamboo charcoal, spandex

Application

- Put the brace on the knee.
- Adjust the strap below the knee cap.

Important Information

- If you have questions regarding application, please contact your local dealer or the prescriber.
- The product is not flameproof, avoid temperature above 120°C / 248°F.
- The function of the product must not be changed or compromised, the product should only be used for its original, intended indication.
- The product is made for possible all day use. For your safety though, see point 5.
- If you are allergic to heat, we advise you to not use the product for more than 3-4 hours at a time.
- Ensure that the support / orthosis is not over-tightened to avoid risk of local pressure and constriction of blood vessels and nerves. Verify the correct positioning and fit of the support / orthosis.
- Do not use the product on more than one patient to avoid skin irritation, formation of eczema or infections.

Washing Instructions

- Close all Velcro.
- Mesh laundry bag is recommended.
- Stretch product while wet and let dry.

Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by the manufacturer with sole responsibility according to appendix VII of the guidelines.

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Material

Polyamid, polyester/bamboo charcoal, spandex

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Otto Bock HealthCare GmbH