


KNEES Keep the knee slightly bent (about 30°) and measure the circumference (cm) around the leg, 10 cm below the patella.



	XS	S	M	L	XL	XXL
6903		34-36	37-39	40-42		
7750, 7751, 7752, 7754, 7763, 7765, 7781, 7782, 7783	31-33	33-35	35-37	37-40	40-43	43-46
7751W, 7753W	30-32	32-34	34-36	36-38	38-40	
7757	27-35		35-43			
7759, 7790		33-35	35-37	37-40	40-43	43-46
7952		26-27	27-28	28-30	30-32	
7953	28-30	30-33	33-36	36-39	39-42	